


**Where to measure:**

The iFIT transfemoral prosthesis will fit limb lengths 17cm to 30+cm measured from adductor tubercle to end of residual limb. Please refer to instruction guide for detailed instructions.



**Sizing Chart:**

Size	Limb Circumferences	
	Distal (cm)*	Proximal (cm)*
Standard	30 – 44 cm	36 – 50 cm
Wide	43 – 60 cm	50 – 65cm

\*These measurements are taken while the patient is wearing a 3mm liner. There may be some overlap due to compressibility of the residual limb soft tissues. Also take into account some patients may have a bulbous or conical shaped limb.

**Transfemoral Prosthesis Images:**

